

Micro-volume UV/Vis Spectrophotometer

UVISDrop

UVISDrop gives highly accurate full spectrum measurements with micro-samples (spectrum range: 200~ 850 nm). Simply pipette 0.5~ 2.5ul sample onto the pedestal, laying down the sampling arm, UVISDrop can automatically complete the measurement within 5 seconds.

Features

1. UV/Vis measurement. Wide application, including nucleic acid (DNA/ RNA), protein, the concentration of other samples and counting of cell numbers.
2. A full spectrum measurement (200~ 850 nm) is available. The results of two light paths can be displayed simultaneously.
3. No need to dilute high-concentration samples.
4. No cuvettes or capillaries or other positioning instruments required.
5. Buttons of "Blank" & "Measurement" on the shell, no need to use mouse.
6. The light intensity is adjustable according to the sample concentrations to obtain stable measurements.
7. Controllable by the computer software with data storage function.

Specification

Standard path length	1 mm
Short path length (for high concentration sample)	0.2 mm
Sample volume	0.5 ~ 2.5ul
Light source	Xenon flash lamp
Detector	3648-element linear silicon CCD array detector
Wavelength range	200~ 850 nm
Wavelength accuracy	1 nm
Wavelength resolution	3 nm (FWHM at Hg 546 nm) ≤ 1.8 nm (FWHM at Hg 253.7 nm)
Absorbance precision	0.002 (1 mm path)
Absorbance accuracy	2% (at 0.76 absorbance at 257 nm)
Measurement range	0.1 ~ 70 Abs.
Measurement accuracy	0.1 ~ 5 Abs. ±0.1 5 ~ 75 Abs. ±2%
Measurement cycle	< 5 seconds
Size	214 x147 x 95 mm
Optical fiber	Quartz
Voltage	12 VDC
Power	6W
Detects nucleic acid up to	2-3500 ng/ ul
Detects protein up to	100 mg/ ml
Weight	1.4 kg



UVS-99